

PDS5A antibody

Catalog Number: 113637

Product name

PDS5A antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDS5A Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human colon tissue. Positive WB detected in HeLa cells, COLO 320 cells, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 130kd

Preparation

This antibody was obtained by immunization of PDS5A recombinant protein (Accession Number: BC009650). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113637(PDS5A antibody) at dilution of 1:400



Immunofluorescent analysis of HeLa cells, using PDS5A antibody Catalog No:113637 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).




Immunohistochemistry of paraffin-embedded

PDS5A antibody

Catalog Number: 113637

human colon tissue slide using Catalog No:113637(PDS5A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113637(PDS5A Antibody) at dilution of 1:50 (under 40x lens)