

# PDXDC1 antibody



Catalog Number: 113641

## Product name

PDXDC1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PDXDC1 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis. Positive IF detected in HepG2. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 87kd

## Preparation

This antibody was obtained by immunization of PDXDC1 recombinant protein (Accession Number: BC033748). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113641(PDXDC1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:113641(PDXDC1 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded

# PDXDC1 antibody

Catalog Number: 113641

---

human hepatocirrhosis using Catalog  
No:113641(PDXDC1 antibody) at dilution of 1:100  
(under 40x lens)



Immunofluorescent analysis of HepG2 cells, using  
PDXDC1 antibody Catalog No:113641 at 1:25  
dilution and Rhodamine-labeled goat anti-rabbit  
IgG (red).