

PDZD3 antibody

Catalog Number: 113644

Product name

PDZD3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDZD3 Rabbit Polyclonal antibody. Positive WB detected in rat kidney tissue, mouse colon tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 54-60 kDa

Preparation

This antibody was obtained by immunization of PDZD3 recombinant protein (Accession Number: NM_024791). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



rat kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113644(PDZD3 Antibody) at dilution of 1:600



IP Result of anti-PDZD3 (IP:Catalog No:113644, 4ug; Detection:Catalog No:113644 1:500) with mouse kidney tissue lysate 4800ug.