

PDCL2 antibody

Catalog Number: 113653

Product name

PDCL2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PDCL2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IHC detected in human kidney. Observed molecular weight by Western-blot: 30-35kd

Preparation

This antibody was obtained by immunization of PDCL2 recombinant protein (Accession Number: BC034431). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113653(PDCL2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113653(PDCL2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113653(PDCL2 antibody) at dilution of 1:50 (under 40x lens)