

# PDE10A antibody

Catalog Number: 113655

## Product name

PDE10A antibody

## Immunogen

[Human PDE10A Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDE10A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 89 kDa

## Preparation

This antibody was obtained by immunization of PDE10A recombinant protein (Accession Number: XM\_047418099). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IP

## Dilutions

Recommended Dilution:  
IP: 1:200-1:1000

## Validations



IP Result of anti-PDE10A (IP:Catalog No:113655, 4ug; Detection:Catalog No:113655 1:300) with mouse brain tissue lysate 5200ug.