

# PDE2A antibody

Catalog Number: 113659

## Product name

PDE2A antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDE2A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, mouse heart tissue. Positive IHC detected in human kidney tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 100-106 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002599). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113659(PDE2A antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113659(PDE2A antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human kidney using Catalog No:113659(PDE2A

# PDE2A antibody

Catalog Number: 113659

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-PDE2A (IP:Catalog No:113659, 4ug; Detection:Catalog No:113659 1:500) with mouse brain tissue lysate 3440ug.