

PDE6C antibody

Catalog Number: 113667

Product name

PDE6C antibody

Specificity

Human; other species not tested.

Antibody description

PDE6C Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 52 kDa; 99 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006204). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:113667(PDE6C antibody) at dilution of 1:200