

# PDE7A antibody



Catalog Number: 113668

## Product name

PDE7A antibody

## Specificity

Human; other species not tested.

## Antibody description

PDE7A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of PDE7A recombinant protein (Accession Number: NM\_002603). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113668(PDE7A Antibody) at dilution of 1:1000