

PDF antibody



Catalog Number: 113673

Product name

Polyclonal

PDF antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

PDF Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 60 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of PDF recombinant protein (Accession Number: bc019912). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113673(PDF Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality