

PDGFRB antibody



Catalog Number: 113676

Product name

PDGFRB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDGFRB Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive FC detected in A549 cells. Observed molecular weight by Western-blot: 160-180 kDa

Preparation

This antibody was obtained by immunization of PDGFRB recombinant protein (Accession Number: NM_002609). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113676(PDGFRB antibody) at dilution of



1:1000

1X10⁶ A549 cells were stained with 0.2ug PDGFRB antibody (Catalog No:113676, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.