

# PDHX antibody



Catalog Number: 113680

## Product name

PDHX antibody

IP: 1:2000-1:20000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDHX Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HEK-293 cells, SKOV-3 cells. Positive IP detected in HEK-293 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of PDHX recombinant protein (Accession Number: NM\_003477). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, IP, WB

## Dilutions


Recommended Dilution:

WB: 1:2000-1:20000


## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113680(PDHX antibody) at dilution of 1:5000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113680(PDHX antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-PDHX (IP:Catalog No:113680,

# PDHX antibody



Catalog Number: 113680

---

4ug; Detection:Catalog No:113680 1:5000) with HEK-293 cells lysate 2000ug.