

# PCDH1 antibody



Catalog Number: 113688

## Product name

Polyclonal

PCDH1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

PCDH1 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, human brain tissue. Observed molecular weight by Western-blot: 115 kDa

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of PCDH1 recombinant protein (Accession Number: NM\_002587). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113688(PCDH1 Antibody) at dilution of 1:300