

# PCDH1 antibody



Catalog Number: 113688

## Product name

PCDH1 antibody

## Immunogen

[Human PCDH1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PCDH1 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, human brain tissue. Observed molecular weight by Western-blot: 115 kDa

## Preparation

This antibody was obtained by immunization of PCDH1 recombinant protein (Accession Number: NM\_002587). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113688(PCDH1 Antibody) at dilution of 1:300