

PCDH11Y antibody



Catalog Number: 113691

Product name

PCDH11Y antibody

Specificity

Human; other species not tested.

Antibody description

PCDH11Y Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 114-147 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_032973).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113691(PCDH11Y antibody) at dilution of 1:300