

# PCDHA2 antibody



Catalog Number: 113693

## Product name

PCDHA2 antibody

## Immunogen

[Human PCDHA2 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PCDHA2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, human heart tissue, Raji cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of PCDHA2 recombinant protein (Accession Number: NM\_031496). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113693(PCDHA2 antibody) at dilution of 1:800



IP Result of anti-PCDHA2 (IP:Catalog No:113693, 3ug; Detection:Catalog No:113693 1:300) with HeLa cells lysate 4650ug.