

# PCDHA6 antibody



Catalog Number: 113696

## Product name

PCDHA6 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PCDHA6 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, K-562 cells. Observed molecular weight by Western-blot: 75kd, 103kd

## Preparation

This antibody was obtained by immunization of PCDHA6 recombinant protein (Accession Number: NM\_018909). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113696(PCDHA6 antibody) at dilution of 1:500