

PCDHA9 antibody

Catalog Number: 113697

Product name

PCDHA9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCDHA9 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 91 kDa

Preparation

This antibody was obtained by immunization of PCDHA9 recombinant protein (Accession Number: NM_031857). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113697(PCDHA9 Antibody) at dilution of 1:500



IP Result of anti-PCDHA9 (IP:Catalog No:113697, 4ug; Detection:Catalog No:113697 1:300) with mouse brain tissue lysate 4000ug.