

PCDHB5 antibody



Catalog Number: 113699

Product name

PCDHB5 antibody

Specificity

human, mouse; other species not tested.

Antibody description

PCDHB5 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells, mouse brain tissue. Observed molecular weight by Western-blot: 70-75kd

Preparation

This antibody was obtained by immunization of PCDHB5 recombinant protein (Accession Number: NM_015669). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113699(PCDHB5 antibody) at dilution of 1:800