

Pepsinogen II antibody



Catalog Number: 113709

Product name

Polyclonal

Pepsinogen II antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

Pepsinogen II Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue. Observed molecular weight by Western-blot: 42kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of Pepsinogen II recombinant protein (Accession Number: NM_002630). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:113709(PGC antibody) at dilution of 1:800