

# PER1 antibody



Catalog Number: 113711

## Product name

PER1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PER1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in A549 cells, human brain tissue, mouse kidney tissue, mouse liver tissue. Observed molecular weight by Western-blot: 50-55kd

## Preparation

This antibody was obtained by immunization of PER1 recombinant protein (Accession Number: BC028207). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:113711(PER1 antibody) at dilution of 1:300



Immunofluorescent analysis of HeLa cells, using PER1 antibody Catalog No:113711 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).