

# PER2 antibody

Catalog Number: 113712

## Product name

PER2 antibody

## Specificity

Human; other species not tested.

## Antibody description

PER2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, K-562 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 136kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_022817). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113712(PER2 antibody) at dilution of 1:100



IP Result of anti-PER2 (IP:Catalog No:113712, 3ug; Detection:Catalog No:113712 1:300) with HEK-293 cells lysate 1000ug.