

# PEX7 antibody



Catalog Number: 113732

## Product name

PEX7 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PEX7 Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of PEX7 recombinant protein (Accession Number: NM\_000288). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:113732(PEX7 antibody) at dilution of 1:800