

PFDN1 antibody



Catalog Number: 113733

Product name

PFDN1 antibody

Specificity

human; other species not tested.

Antibody description

PFDN1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human colon tissue. Positive IHC detected in human colon cancer. Observed molecular weight by Western-blot: 14 kDa

Preparation

This antibody was obtained by immunization of PFDN1 recombinant protein (Accession Number: NM_002622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113733(PFDN1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113733(PFDN1 antibody) at dilution of 1:50 (under 10x lens)