

# PFDN1 antibody



Catalog Number: 113733

## Product name

PFDN1 antibody

## Specificity

human; other species not tested.

## Antibody description

PFDN1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human colon tissue. Positive IHC detected in human colon cancer. Observed molecular weight by Western-blot: 14 kDa

## Preparation

This antibody was obtained by immunization of PFDN1 recombinant protein (Accession Number: NM\_002622). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113733(PFDN1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113733(PFDN1 antibody) at dilution of 1:50 (under 10x lens)