

PFDN2 antibody

Catalog Number: 113734

Product name

PFDN2 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

PFDN2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells. Positive IHC detected in human lung cancer tissue. Positive IF detected in L02 cells. Observed molecular weight by Western-blot: 17-20 kDa




Preparation

This antibody was obtained by immunization of PFDN2 recombinant protein (Accession Number: NM_012394). Purification method: Antigen affinity purified.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113734(PFDN2 antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113734(PFDN2 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF



Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Immunofluorescent analysis of L02 cells using

PFDN2 antibody



Catalog Number: 113734

Catalog No:113734(PFDN2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)