

PFDN5 antibody

Catalog Number: 113736

Product name

PFDN5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PFDN5 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse pancreas tissue, PC-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of PFDN5 recombinant protein (Accession Number: BC003373). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113736(PFDN5 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:113736(PFDN5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)