

PCGF6 antibody



Catalog Number: 113737

Product name

PCGF6 antibody

Specificity

Human; other species not tested.

Antibody description

PCGF6 Rabbit Polyclonal antibody. Positive WB detected in PC-13 cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of PCGF6 recombinant protein (Accession Number: NM_001011663). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:113737(PCGF6 Antibody) at dilution of 1:1500