

# PCIF1 antibody

Catalog Number: 113738

## Product name

PCIF1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PCIF1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human stomach tissue. Positive WB detected in Jurkat cells, HepG2 cells. Observed molecular weight by Western-blot: 81 kDa

## Preparation

This antibody was obtained by immunization of PCIF1 recombinant protein (Accession Number: XM\_011528981). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113738(PCIF1 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:113738(PCIF1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog

# PCIF1 antibody



Catalog Number: 113738

---

No:113738(PCIF1 Antibody) at dilution of 1:50 (under 40x lens)