

PCIF1 antibody



Catalog Number: 113738

Product name

PCIF1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PCIF1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human stomach tissue. Positive WB detected in Jurkat cells, HepG2 cells. Observed molecular weight by Western-blot: 81 kDa

Preparation

This antibody was obtained by immunization of PCIF1 recombinant protein (Accession Number: XM_011528981). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113738(PCIF1 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:113738(PCIF1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog

PCIF1 antibody



Catalog Number: 113738

No:113738(PCIF1 Antibody) at dilution of 1:50 (under 40x lens)