

PCMT1 antibody



Catalog Number: 113742

Product name

PCMT1 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCMT1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human testis tissue, Raji cells. Positive IHC detected in human pancreas tissue. Observed molecular weight by Western-blot: 25kd

Validations




Preparation

This antibody was obtained by immunization of PCMT1 recombinant protein (Accession Number: NM_001360452). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113742(PCMT1 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human pancreas using Catalog No:113742(PCMT1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:113742(PCMT1 antibody) at dilution of 1:50 (under 40x lens)