

PCNXL2 antibody



Catalog Number: 113746

Product name

PCNXL2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PCNXL2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 237-260kd

Preparation

This antibody was obtained by immunization of PCNXL2 recombinant protein (Accession Number: BC008300). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113746(PCNXL2 antibody) at dilution of 1:500