

# PCSK7 antibody

Catalog Number: 113753

## Product name

PCSK7 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PCSK7 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of PCSK7 recombinant protein (Accession Number: NM\_004716). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113753(PCSK7 antibody) at dilution of 1:50



(under 10x lens)

mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113753(PCSK7 Antibody) at dilution of 1:300