

# PCTAIRE3 antibody



Catalog Number: 113756

## Product name

PCTAIRE3 antibody

Polyclonal

## Ig Type

## Specificity

Human, Mouse, Rat; other species not tested.

Rabbit IgG

## Applications

## Antibody description

PCTAIRE3 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 54 kDa

ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of PCTAIRE3 recombinant protein (Accession Number: NM\_212502). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113756(PCTK3 antibody) at dilution of 1:500

## Clonality