

PCYOX1 antibody



Catalog Number: 113758

Product name

PCYOX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCYOX1 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue.

Preparation

This antibody was obtained by immunization of PCYOX1 recombinant protein (Accession Number: NM_016297). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon using Catalog No:113758(PCYOX1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:113758(PCYOX1 antibody) at dilution of 1:50 (under 40x lens)