

PCYT1B antibody

Catalog Number: 113759

Product name

PCYT1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCYT1B Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue, human brain tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of PCYT1B recombinant protein (Accession Number: NM_001163264). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113759(PCYT1B antibody) at dilution of 1:500



IP Result of anti-PCYT1B (IP:Catalog No:113759, 4ug; Detection:Catalog No:113759 1:500) with human placenta tissue lysate 4000ug.