

# PFTK1 antibody



Catalog Number: 113774

## Product name

PFTK1 antibody

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PFTK1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse brain tissue, SH-SY5Y cells. Positive IHC detected in human kidney tissue, mouse brain tissue. Observed molecular weight by Western-blot: 48-53kd

## Preparation

This antibody was obtained by immunization of PFTK1 recombinant protein (Accession Number: NM\_012395). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113774(PFTK1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113774(PFTK1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113774(PFTK1 antibody) at dilution of 1:100 (under 40x lens)