

PGK2 antibody



Catalog Number: 113785

Product name

Polyclonal

PGK2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

PGK2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HeLa cells, HepG2 cells, human heart tissue, human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 45 kDa

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Preparation

Validations

This antibody was obtained by immunization of PGK2 recombinant protein (Accession Number: NM_138733). Purification method: Antigen affinity purified.


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113785(PGK2 antibody) at dilution of 1:1200