

PGLYRP3 antibody



Catalog Number: 113787

Product name

Polyclonal

PGLYRP3 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

PGLYRP3 Rabbit Polyclonal antibody. Positive WB detected in Apoptosised HeLa cells, mouse brain tissue, mouse colon tissue. Observed molecular weight by Western-blot: 38kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of PGLYRP3 recombinant protein (Accession Number: XM_011509119). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Apoptosised HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113787(PGLYRP3 antibody) at dilution of 1:500