

PGT antibody



Catalog Number: 113796

Product name

Polyclonal

PGT antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

PGT Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, mouse heart tissue, mouse liver tissue. Observed molecular weight by Western-blot: 65 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of PGT recombinant protein (Accession Number: NM_005630). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113796(PGT antibody) at dilution of 1:1000

Clonality