

# PHF17 antibody



Catalog Number: 113809

## Product name

PHF17 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PHF17 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 58-61kd, 94-96kd

## Preparation

This antibody was obtained by immunization of PHF17 recombinant protein (Accession Number: NM\_024900). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113809(PHF17 antibody) at dilution of 1:500