

PHF21A antibody



Catalog Number: 113813

Product name

Polyclonal

PHF21A antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

PHF21A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 75 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of PHF21A recombinant protein (Accession Number: NM_001101802). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113813(PHF21A antibody) at dilution of 1:500