

# PHF6 antibody



Catalog Number: 113815

## Product name

PHF6 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PHF6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Positive IHC detected in human ovary tumor tissue, human cervical cancer tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 41 kDa

## Preparation

This antibody was obtained by immunization of PHF6 recombinant protein (Accession Number: NM\_032335). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113815(PHF6 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:113815(PHF6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded


# PHF6 antibody



Catalog Number: 113815

---

human ovary tumor using Catalog  
No:113815(PHF6 antibody) at dilution of 1:50  
(under 40x lens)



Immunofluorescent analysis of (10%  
Formaldehyde) fixed HeLa cells using Catalog  
No:113815(PHF6 Antibody) at dilution of 1:50 and  
Alexa Fluor 488-conjugated AffiniPure Goat Anti-  
Rabbit IgG(H+L)