

PHOX2A antibody



Catalog Number: 113831

Product name

PHOX2A antibody

Specificity

Human; other species not tested.

Antibody description

PHOX2A Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in SH-SY5Y cells, HeLa cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of PHOX2A recombinant protein (Accession Number: NM_005169). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113831(PHOX2A Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed SH-SY5Y cells using Catalog No:113831(PHOX2A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)