

PKA C-alpha specific antibody



Catalog Number: 113844

Product name

PKA C-alpha specific antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

PKA C-alpha specific Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human testis tissue, human liver tissue. Positive IHC detected in human colon cancer tissue, human testis tissue. Observed molecular weight by Western-blot: 38-43 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002730). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113844(PRKACA antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human colon cancer slide using Catalog No:113844(PRKACA Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human colon cancer slide using Catalog

PKA C-alpha speicifc antibody



Catalog Number: 113844

No:113844(PRKACA Antibody) at dilution of 1:50

IP Result of anti-PRKACA (IP:Catalog No:113844, 4ug; Detection:Catalog No:113844 1:500) with mouse brain tissue lysate 2640ug.

