

PKD2L1 antibody



Catalog Number: 113854

Product name

PKD2L1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PKD2L1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 92kd

Preparation

This antibody was obtained by immunization of PKD2L1 recombinant protein (Accession Number: NM_016112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113854(PKD2L1 antibody) at dilution of 1:500