

# PHYHD1 antibody



Catalog Number: 113865

## Product name

Polyclonal

PHYHD1 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

PHYHD1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue, L02 cells. Observed molecular weight by Western-blot: 30-32kd

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of PHYHD1 recombinant protein (Accession Number: NM\_174933). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113865(PHYHD1 antibody) at dilution of 1:800