

PI16 antibody



Catalog Number: 113868

Product name

Polyclonal

PI16 antibody

Ig Type

Specificity

Rabbit IgG

human; other species not tested.

Applications

Antibody description

ELISA, WB

PI16 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, COLO 320 cells, SKOV-3 cells. Observed molecular weight by Western-blot: 49-60kd,30-35kd

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of PI16 recombinant protein (Accession Number: BC022399). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:113868(PI16 antibody) at dilution of 1:200

Clonality