

# PI3K p55(gamma) antibody

Catalog Number: 113874

## Product name

PI3K p55(gamma) antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PI3K p55(gamma) Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 110 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003629). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113874(PIK3R3 antibody) at dilution of 1:500



IP Result of anti-PIK3R3 (IP:Catalog No:113874, 5ug; Detection:Catalog No:113874 1:800) with mouse testis tissue lysate 4800ug.