

PIAS1 antibody



Catalog Number: 113880

Product name

PIAS1 antibody

Specificity

human; other species not tested.

Antibody description

PIAS1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HepG2 cells.

Preparation

This antibody was obtained by immunization of PIAS1 recombinant protein (Accession Number: NM_016166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113880(PIAS1 antibody) at dilution of 1:1000