

PIAS4 antibody



Catalog Number: 113884

Product name

PIAS4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PIAS4 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HL-60 cells, mouse ovary tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 57-70kd

Preparation

This antibody was obtained by immunization of PIAS4 recombinant protein (Accession Number: NM_015897). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113884(PIAS4 antibody) at dilution of 1:300